



miRBase precursor : ath-MIR156c
 Total read count : 24
 ath-miR156a read count : 4
 ath-miR156b read count : 4
 ath-miR156c read count : 4
 ath-miR156d read count : 4
 ath-miR156e read count : 4
 ath-miR156f read count : 4
 remaining reads : 0



- ath-miR156e
- ath-miR156a
- ath-miR156d
- ath-miR156f
- ath-miR156b
- ath-miR156c

5'	cgcauagaaacugacagaagagagugagcac	acaaaggcacuuugcauguucgaugcauuugcuucucucugcgugcucacugcucuaucugucagauucggcu	-3'	exp	reads	mm	sample
gacagaagagagugagcac.....				4	0	seq