


miRBase precursor : ath-MIRf10256-akr  
 Total read count : 24  
 ath-miR156a read count : 4  
 ath-miR156b read count : 4  
 ath-miR156c read count : 4  
 ath-miR156d read count : 4  
 ath-miR156e read count : 4  
 ath-miR156f read count : 4  
 ath-miRf10256-akr read count 0  
 ath-miRf10257-akr read count 0  
 remaining reads : 0



ath-miR156d  
 ath-miR156b  
 ath-miR156a  
 ath-miR156e  
 ath-miR156f  
 ath-miR156c  
 ath-miRf10257-akr  
 ath-miRf10256-akr

5' -	ggugacagaagagagugagcacacacagguggcuuucugcauuugaagguuccaugcuugaagcuaugugucacacucucuauccgucacc	-3'	exp			
	(((((.....(((((.....)))))))))))))))))))))))))))))	reads	mm	sample		
...	gacagaagagagugagcac.....	4	0	seq		